



Section 1: Clinical Modality Handpieces					
Number	Item Image Context	Name	QTY	Intended Clinical Function & Use Case	Note
1		Premium Handheld Laser Handpiece	1	Delivers high-power laser therapy (Class 4) for deep tissue photobiomodulation, accelerated cellular repair, and deep-seated pain relief.	\$375/pcs
2		Magnetic Therapy Handpiece	1	The core EMTT (Physio Magneto) applicator. Emits high-oscillating magnetic fields to stimulate cell membrane permeability and accelerate musculoskeletal regeneration.	\$450/pcs
3		Cold Laser Cluster Handpiece	1	Designed for Low-Level Laser Therapy (LLLT) or Cold Laser treatments. Targets superficial tissue layers, promotes wound healing, and reduces localized inflammation over broad body areas.	\$225/pcs
4		Reusable Electrotherapy Electrode Pad (with Integrated Lead Wire)	1	Must be used in conjunction with a biocompatible hydrogel self-adhesive interface pad to ensure safe contact and flawless current transfer to the patient's skin.	\$10/pcs

Section 2: Specialized Laser Treatment Tips & Optical Safeties					
Number	Item Image Context	Name	QTY	Intended Clinical Function & Use Case	Note
1		High-Power Laser Therapy Tips	1	Interchangeable spacers that thread onto the laser handpiece to adjust beam spot size, divergence, and treatment distance based on the target area's anatomy.	D10: \$50/pcs D30: \$100/pcs D60: \$90/pcs
2		Quartz Crystal Treatment Head	1	A specialized crystal lens head providing a highly focused, convergent laser beam profile for concentrated point treatments and deep trigger points.	\$88.50/pcs
3		Optical Port Safety Termination Cap	1	Mandatory Safety Device. Must be screwed onto the machine's primary laser port whenever the handpiece is disconnected to prevent accidental beam emission and dust ingress.	
4		Patient Laser Safety Eye Shield	1	Complete optical block goggles designed to protect the patient's eyes from scatter radiation during facial or cervical laser procedures.	
5		Clinician Laser Safety Glasses		Advanced filter glasses worn by the medical practitioner to block harmful laser wavelengths while maintaining clear ambient room visibility.	\$225/pcs

Section 3: Structural Hardware, Brackets & Workstation					
Number	Item Image Context	Name	QTY	Intended Clinical Function & Use Case	Note
1		Primary Touchscreen Display Console	1	The core control screen assembly for monitoring system status and choosing software protocols.	\$525/pcs
2		Heavy-Duty Articulating Support Arm	1	Provides structural support and flexible positioning for the EMTT loop or laser handpieces during passive static therapies.	\$330/pcs
3		Laser Handpiece Cradle Mount	1	Secure click-in holster mounted on the machine frame for stable docking of the laser handpiece.	
4		Anti-Slip Structural Silicone Shim	1	Secure click-in holster mounted on the machine frame for stable docking of the laser handpiece.	

Section 4: Structural Hardware & Assembly Toolset					
Number	Item Image Context	Name	QTY	Intended Clinical Function & Use Case	Note
1		6-inch Adjustable Spanner Wrench	1	Multi-purpose adjustable mechanic tool for custom hardware tightening and fluid line compression fittings.	
2		L-Shape Allen Wrench Set	1	Standard right-angle hex keys for deep-	

3		T-Handle Hex Wrench (H3 / H3T)	1	High-torque T-handle hex wrench specifically optimized for quick and ergonomic adjustment of frame pivot points.	
4		Precision Phillips Screwdriver	1	Standard cross-head screwdriver for mounting electronics housing and secure panel screws.	

Section 5: Water Management System & Fluid Calibration					
Number	Item Image Context	Name	QTY	Intended Clinical Function & Use Case	Note
1		Fluid Refill Tube with Funnel	1	Connects directly to the machine's "Inlet" port to cleanly add purified/distilled water to the system reservoir.	
2		Fluid Drainage Tube	1	Attaches securely to the "Overflow/Outlet" port for safe routine draining or cycling of fluid during maintenance.	
3		Quick-Connect Fluid male shipping plug	1	Upon receiving the machine, this male shipping plug must be discarded. Otherwise, an air pressure imbalance will occur inside the water tank, which may cause the tank to burst.	

The replacement price for the circuit board is \$120. Please note that the shipping fee for a small package currently ranges from \$150 to \$180.